



DARBY DENTAL SUPPLY, LLC

Temporary Crown and Bridge 10:1 Material

Darby Temporary Crown & Bridge 10:1 Material is an enhanced crown and bridge material that offers improved physical properties. This material provides predictable temporization for a range of indirect restorations including single crowns to extended bridges. The rapid automix delivery system assures a perfect chemistry every time.

Features and Benefits:

- 4 Minute Setting Time
- Easy to Trim
- Excellent Physical Properties
- Stain Resistant
- Polishes Easily
- Easy to Extrude
- Patient Friendly

50 ml HP Cartridge,

10 HP Mixing Tips

9430100 - A1

9430102 - A2

9430104 - A3

9430106 - A3.5


9430108 - B1

9430110 - Bleach

10:1 Mixing Tips

25/Bag

9515144



Darby Dental Supply, LLC • 800-645-2310 • www.darby.com