



DARBY DENTAL SUPPLY, LLC

dry tips[®]

ULTIMATE SALIVA CONTROL



IDEAL FOR:

IMPRESSIONS

OCCLUSION
GRINDING

ORTHODONTICS

PAINT-ON
TECHNIQUES

COMPOSITES

CEMENTATION

FISSURE SEALING

FLUORIDE
APPLICATIONS

dry
tips®



SALIVA ABSORPTION CAPACITY 30 TIMES ITS WEIGHT

Far outlasts cotton rolls and other absorbent pads,
and doesn't lose moisture when handled.



- Effective replacement for rubber dams
- Complete coverage of parotid orifice and entire buccal mucosa with polyethylene film for cutting protection.
- Facilitates work far back in oral cavity.
- No loose fibers or other residues left in mouth.
- Flexible; accommodates to cheek movements.
- Comfortably stays in place, yet is easily removed with water spray.

Dry Tips Ordering Information

Item Number	Size	Quantity	Color
9532461 (291543)	Small	50	White
9532463 (291542)	Large	50	White

*Order today through your preferred dental dealer. For questions, call Microbrush International at 1 (262) 375-4011.

microbrush.com

Convenient. Consistent. Control.™

To Reorder Call 800.645.2310 or Visit darby.com



Microbrush