



DARBY DENTAL SUPPLY, LLC



INTRODUCING the **ProphyPal™** HYGIENE HANDPIECE

What's your ProphyPal personality?

Prophy
Pink



Classic
Silver



Wild
Thing
Green



1/Pkg.

8281232 - Prophy Pink

8281230 - Classic Silver

8281234 - Wild Thing Green

3/Pkg.

8281242 - Variety

8281238 - Prophy Pink

8281236 - Classic Silver

8281240 - Wild Thing Green

The ProphyPal comes in 3 great colors to color-code your practice for sterilization cycles, or to add some personality to the office!

Denticator®
SINCE 1931
www.denticator.com

Hygienists Spoke, We Listened

360° Swivel and Grip Nosecone - allows for maximum control and easy transfer between buccal and lingual surfaces


Sure-fit™ Nosecone - extended nosecone provides excellent stability for all disposable prophy angles

Power and Torque - 0 – 5,000 RPM motor operates at optimal speed and torque to provide tough polishing on heavy stain without stalling

Ergonomic - the lightweight and balanced construction provides excellent ergonomic positions for all hand sizes

Durable - built in the USA from quality materials, the ProphyPal comes with a **2 year comprehensive warranty**

Whisper Quiet - air powered operation reduces patient and operator stress



Additional Product Features

- ◇ Designed specifically for disposable prophy angles
- ◇ Autoclavable up to 275°F/135°C
- ◇ 4-hole standard tubing connection
- ◇ Operates using 40-45 psi
- ◇ Proudly manufactured in the U.S.A.

To Reorder Call 800.645.2310
or Visit darby.com

