

MIDWEST®

SKILLFUL HANDLING. POWERFUL DRIVE.

Introducing MIDWEST® *E* – The electric handpiece system that provides a necessary balance of power and control for all your high-speed, low-speed, and endodontic procedures

NEW

MIDWEST® *E*
Electric Handpiece System

NEW

MIDWEST[®] E

Electric Handpiece System



THE NECESSARY BALANCE OF POWER AND CONTROL

E Power—lightweight autoclavable motor weighs only 59 g—among the lightest offered

E Control—touchscreen user interface, small handpiece designed for easy maneuvering, and the ability to make a variety of adjustments on a single screen

E Attachments offered in high speed, low speed and endodontic configurations

E Backed by the Legendary Midwest quality and support



ULTIMATE CONTROL

- Touch screen control offers one-screen review of operating settings
- Endo functionality includes user defined torque limits and 3 auto-reverse options

ATTACHMENTS

SYSTEM COMPONENTS

SKU #	Vendor #	Item
8640358	876000	Midwest [®] Electric Handpiece System Includes: CONTROL UNIT, MOTOR, POWER SUPPLY AND HOSE
8640360	875005	Midwest [®] E High Speed Attachment 1:5 HS FO CA (1:5L)
8640362	875010	Midwest [®] E Attachment 1:1 LS FO CA (1:1L)
8640364	875015	Midwest [®] E Attachment 1:1 LS CA (1:1)
8640366	875020	Midwest [®] E Speed Reducing Attachment 5:1 LS FO CA (5:1L)
8640368	875025	Midwest [®] E Speed Reducing Attachment 5:1 LS CA (5:1)
8640370	875030	Midwest [®] E Straight Attachment 1:1 Straight (1:1S)
8640372	875040	Midwest [®] E Electric Control Unit Only
8640374	875045	Midwest [®] E Electric Motor Only
8640376	875050	Midwest [®] E Electric Motor Hose Only
8640378	875080	Midwest [®] E Remote Mounting Kit



DARBY DENTAL SUPPLY, LLC

To Reorder Call
800.645.2310
or Visit darby.com

